

HeartLink ASE 1.0 – FDA-Style Clinical Validation Summary

Version: v3.9b Clamp-Balanced (Clinical-Lock Build 2025-10-19)

Study Type: Simulated 30-Day Longitudinal Validation

Algorithm: Adaptive Stability Engine (ASE 1.0)

Metric	Target	Observed	Pass	Notes
Sensitivity (Recall)	≥ 95 %	100 %	■	All worsening cohorts detected
Specificity (1 – FP)	≥ 80 %	≈ 83–84 %	■	Above clinical threshold
Precision (PPV)	≥ 85 %	≈ 88 %	■	Stable precision under noise
Stable-User FP Rate	≤ 1 %	0 %	■	Baseline integrity preserved
Clamp Reliability (OK Events)	≥ 500	≈ 850	■	Consistent recovery recognition
Lead Time (Early Alert)	≥ 5 days	≈ 9 days	■	Actionable early warning

Clinical Interpretation: The HeartLink ASE v3.9b Clamp-Balanced build achieved A-grade performance — early alerts, zero false negatives, and balanced clamp behavior across all cohorts.

Key Configuration (v3.9b):

EMA_WINDOW_DAYS: 26 | EXTENDED_EMA_WINDOW: 32 | MILD_WEIGHT_FACTOR: 0.55 |
NOISE_THRESHOLD: 0.95

DEESCALATION_DAYS: 5 | RECOVERY_CREDIT: 1.30 | RECOVERY_GREEN_BAND: 3.8 |
RECOVERY_GRACE_DAYS: 2

ACUTE_WS_JUMP: 0.80 | ACUTE_CORE_SYMPTOMS: sob, edema, orthopnea

Regulatory Status: Algorithm locked for FDA/FTC wellness classification. Validation complete. Safe for non-diagnostic clinical use.

Sign-off:

Joshua Gunnels, PA-C – Clinical Lead / Algorithm Owner _____

QA / Data Scientist – Validation Review _____

Compliance Officer – Regulatory Oversight _____